

# PRODUCT COMPARISON CHART



Your Vision, Our Future

<b>EBUS Bronchoscopes</b>	<b>Universal Platform</b>	<b>NBI</b>	<b>HD</b>	<b>Wide Angle (170)</b>	<b>Close Focus (2mm)</b>
BF-UC160F-OL8	x				
<b>Fiber-optic Bronchoscopes</b>					
BF-XP60	•				
BF-MP60	•				
BF-P60	•				
BF-1T60	•				
BF-3C40	•				
BF-XT40	•				
<b>Video Bronchoscopes</b>					
BF-Q180	x	x			
BF-Q180-AC	x	x			
BF-P180	x	x			
BF-1T180	x	x			
LTF-160	x				
BF-XP160F	x				
BF-MP160F	x				
BF-3C160	x				
BF-XT160	x				
<b>Capsule Endoscopes</b>					
Endo Capsule	x				
<b>Choledochoscopes</b>					
CHF-BP30	•				
<b>Colonoscopes:</b>					
CF-H180AL/I	x	x	x	x	x
CF-Q180AL/I	x	x		x	
PCF-Q180AL/I	x	x			
CF-Q160ZL/I	x	x			
PCF-H180AL/I	x	x	x		x
CF-2T160L/I	x				
<b>Duodenoscopes</b>					
TJF-160VF	x				
<b>Echoendoscopes</b>					
GF-UE160-AL5	x				
GF-UM160	x				
GF-UC140P-AL5 / GF-UCT140-AL5	x				
GF-UC160P-OL5 / GF-UCT160-OL5	x				

• **Fiberscopes - can connect to the CV-180 with coupler adaptor and camera head**

# PRODUCT COMPARISON CHART



Your Vision, Our Future

	<b>Universal Platform</b>	<b>NBI</b>	<b>HD</b>	<b>Wide Angle (170)</b>	<b>Close Focus (2mm)</b>
<b>Enteroscopes</b>					
SIF-Q180	x	x			
<b>Esophagosopes</b>					
PEF-V	x				
<b>Gastrosopes</b>					
GIF-H180	x	x	x		x
GIF-Q180	x	x			
GIF-N180	x	x			
GIF-XP160	x				
GIF-1TQ160	x				
GIF-2T160	x				
GIF-XTQ160	x				
GIF-Q160Z	x	x			
<b>Laryngoscopes</b>					
LF-GP	•				
LF-TP	•				
LF-DP	•				
<b>Rhinolaryngoscopes</b>					
ENF-V2	x	x			
ENF-VQ	x	x			
ENF-VT2	x	x			x
ENF-GP	•				
ENF-P4	•				
ENF-XP	•				
<b>Sigmoidoscopes</b>					
PCF-S	x				
CF-Q160S	x				
OSF-V60					
OSF-4					
<ul style="list-style-type: none"> <li>• <b>Fiberscopes - can connect to the CV-180 with coupler adaptor and camera head</b></li> </ul>					